

# SINGLE USE SUPPORT.

PIONEERING BIOPHARMA



# RoSS.SHIP MULTI-USE

Smart Cold Chain Shipping Container.



[www.susupport.com](http://www.susupport.com)

# RoSS.SHIP Multi-Use


Smart Cold Chain Shipping Container.


**RoSS.SHIP Multi-Use** is a device for storing and shipping RoSS® Shells (Robust Storage and Shipping) containing frozen liquids. The platform is made for carrying large quantities of dry ice, making it the ideal prerequisite for this application.

Article Nr.	Net weight	Material	Cooling Tech.	RoSS Types	RoSS Capacity	Outer dimensions (l x d x h)	Temp. Range	Temp. Duration
SHIP_10xx	200 kg	PUR/PA6	passive - dry ice	all RoSS possible up to 20 liters*	depending on size and RoSS type	1200 x 1200 x 970mm	to -70°C	min. 144 h



Shipment Summary Report RoSS.SHIP Multi-Use



 **Single Use Support assists you from beginning to end, offering highly qualified support all along the way. Contact us for further details!**