


SHIPMENT SERVICE OVERVIEW



SKIP SHIPMENT TROUBLES

Single Use Support takes care of your full shipment.

- Our service package contains all services from pickup to safe delivery.
- Our staff guarantees reliable shipment operations.



CONTACT US AND [CLICK HERE](#)



LEARN MORE ABOUT RoSS.SHIP [HERE](#)

SINGLE USE SUPPORT = ONE STOP SHOP

RoSS.SHIP Single-Use



RoSS.SHIP Multi-Use



ULTRA SMART TRACKING

SMART SURVEILLANCE OF YOUR DRUG SUBSTANCE

Tracking Details	
Temperature	Ranging from -200°C to +80°C
Geo-Position	24/7. GPS, GLONASS, GALILEO, WIFI, GSM)
G-Force	Accelerometer adjustable threshold from 2 g to 12 g
Air Pressure	Range from 300 to 1.100 hPa
Humidity	Definable threshold
Dimensions	129 x 129 x 35 mm



SKIP QUALIFICATION TROUBLES

- **Certified:** 21 CFR Part 11 tracking of your drug substance on the way
- **Conformity:** Calibrated to GDP guideline
- **Cloud-based** storage of tracked data
- **Automated** alarms, reports & archiving



CONTACT US AND CLICK [HERE](#)



LEARN MORE ABOUT RoSS.SHIP [HERE](#)

SINGLE USE SUPPORT = ONE STOP SHOP

Smart Tracking for SHIP

